

# Connex®VM Competency Test

1. True or False.  
Welch Allyn Connex®VM is a data management system that allows a user to collect and review patient data and communicate that data to an information system.
2. The steps to log on to Connex®VM are:
  - a. enter Password, click Login, then enter User Name
  - b. click Login, enter User Name, and Password
  - c. enter User Name and Password, then click Log on
  - d. display Home Page, enter Password, enter User Name
3. True or False.  
After logging on, the Home page appears with a list of patients and their most recent readings with date and time. The list of patients is determined by the unit(s) selected in the Show list.
4. To capture vitals from an attached device, log on, select a patient and click on ▼ next to patient name. Click on Capture vital signs, then Start to obtain vitals, and then Save which will return the user to:
  - a. Log off
  - b. Home page
  - c. Patient Record
  - d. My locations
5. The Capture vital signs screen lets the user document all data on:
  - a. one screen
  - b. one book
  - c. one Patient Note
  - d. My Reminders
6. The purpose of the Refresh button on the Home page is to:
  - a. Sort the patients alphabetically
  - b. Close the Home page
  - c. Update the display to reflect recent changes
  - d. Sort the patient list by room number
7. True or False.  
By clicking on a column heading on the Home page, the user can change the order of the information displayed.

8. To manually enter vitals on the **Capture vital signs** screen that would usually be generated by a vital signs device, you must click on \_\_\_\_\_, then type the values into Blood Pressure, Pulse Rate, Temperature, or SpO<sub>2</sub> fields.
  - a. Enter manually
  - b. Patient Notes field
  - c. Save button
  - d. My Patients
  
9. In order for a clinician to assign or unassign patients to themselves, they must click on:
  - a. Done
  - b. Assign/Unassign
  - c. My Patients
  - d. Save
  
10. True or False.  
To view saved vital signs, from the Home page click on ▼ next to a patient name and select **View patient record**. You can also obtain a printout of the current view by clicking the printer icon.
  
11. To search for a patient not listed in the Patient list, the user would:
  - a. click Home, click Go, type in Last name, First name, or Patient ID
  - b. click Search tab or click on View on toolbar, then click Search. Click on Patients, enter your search criteria, then click on Search.
  - c. click Home, Address Book icon, Done
  - d. none of the above
  
12. .When you are finished charting and you want to exit the application and logout, click on:
  - a. Login
  - b. My Patients
  - c. Log off
  - d. Help

# Connex®VM Competency Test

## Answer Sheet

Name: \_\_\_\_\_

Unit: \_\_\_\_\_

Date: \_\_\_\_\_

1. \_\_\_\_\_

2. \_\_\_\_\_

3. \_\_\_\_\_

4. \_\_\_\_\_

5. \_\_\_\_\_

6. \_\_\_\_\_

7. \_\_\_\_\_

8. \_\_\_\_\_

9. \_\_\_\_\_

10. \_\_\_\_\_

11. \_\_\_\_\_

12. \_\_\_\_\_

# Connex®VM Competency Test

## Answer Key

1. True

2. c

3. True

4. b

5. a

6. c

7. True

8. a

9. b

10. True

11. b

12. c