



Lamps

LED Replacement Lamps	208
Replacement Lamps	210

Batteries

Replacement Batteries	216
-----------------------------	-----

Blood Pressure

Cuff Reference Guide	217
Replacement Connectors	220

Tubing and Accessories

Replacement Tubes and Accessories	222
---	-----

LED LAMPS

REPLACEMENT PARTS



LED Replacement Lamps for Physical Assessment Instruments

While LED technology has been on the market for mainstream use for years, our engineers have worked to optimize its benefits for physical assessment applications by combining color selection and patented circuitry to preserve the feel of halogen while offering the longevity and efficient performance of LED technology.

- Long life backed by a 7-year warranty ensuring peace of mind and virtually eliminating inconvenient interruptions for lamp replacement
- SureColor™ LED technology enabling reduction of intensity to alleviate glare or enhance patient comfort while preserving the color and clarity of the view
- High energy efficiency resulting in 3X the examination time with each charge of a Welch Allyn rechargeable handle
- A simple solution presented in four lamp upgrade kits compatible with our most popular physical assessment instruments



03100
LED



06500
LED



03800
LED

03100-LED Traditional Diagnostic Oscopes

03100-LED Box of 10

06500-LED MacroView™ Diagnostic Oscopes

06500-LED Box of 10

04900-LED Traditional Ophthalmoscopes

04900-LED Box of 10

03800-LED PanOptic™ Ophthalmoscopes

03800-LED Box of 10

A Simple, Flexible Solution

Four lamp upgrade kits cover a broad range of physical assessment instruments you use every day.

Lamp bases are date-stamped for easy warranty verification. Kits include a lamp removal tool for easier replacement.



Dimming Cycle

LED



Consistent color output through dimming cycle





@100%

@25%


Halogen





































LED Lamps

LED Lamp Part #	Device Part #	Description
03100-LED 	Traditional Diagnostic Oscopes	
	20000	Otoscope/Throat Illuminator
	20200	Pneumatic Otoscope
	21700	Operating Otoscope
	25020	Diagnostic Otoscope (w/ Specula)
06500-LED 	MacroView™ Diagnostic Oscopes	
	23810	MacroView™ Otoscope
	23820	MacroView™ Otoscope/Throat Illuminator
	23920	Digital MacroView Otoscope
04900-LED 	Traditional Ophthalmoscopes	
	11720	3.5 V Coaxial Ophthalmoscope
	11730	AutoStep® Coaxial Ophthalmoscope
	11735	Prestige™ Coaxial-Plus Ophthalmoscope
03800-LED 	PanOptic™ Ophthalmoscopes	
	11810	PanOptic™ Ophthalmoscope
	11820	PanOptic™ Ophthalmoscope/Blue Filter/Lens













































LED Color Temperature Spectrum





















Replacement Lamps

Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
00200-U		<p>Illuminators, Transilluminators</p>  41000  26030 <p>Anoscopes</p>  38800  380 Series  385  395	<p>Otoscope</p>  235 <p>PenLite</p>  77900  26000 <p>Illuminators</p>  20100  21600 <p>Transilluminators</p>  42700  40510
00300-U		 60835	
00600-U			<p>Sigmoidoscope</p>  30800
00900-U			<p>Ophthalmoscopes</p>  10600  121  110  17600  17610
01200-U		<p>Binocular Indirect Ophthalmoscope</p>  125 Series	
01300-U			<p>Tongue Blade Holder</p>  28000 <p>PenLites</p>  76000  77700
01400-U			<p>PenLite</p>  77800
02000-U			<p>Headlights</p>  17200












































Replacement Lamps

Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
00200-U			Headlight 46003 
02600-U			Headlight 49003 
03000-U		Ophthalmoscope StrabismoScope EpiScope  11710  12400  47300	Ophthalmoscopes Retinoscope  11600  11605  11610  18000
03100-U		Otosopes  25020  21700  20200  20202 3.5 V Handle Adapter Tongue Blade Holder  73500  28100 Illuminators/Transilluminators  27000  27050  26530  41100  43300	Otosopes  25000  25200  20000
03300-U			Ophthalmoscopes  11511  11500
03400-U			Otosopes  24011  24020  24000  21110, 21111 Otoscope 2.5 V Handle Adapter Illuminators Transilluminator  24031  73550  27200  27250  41110
03600-U		Disposable Laryngoscope Illuminator 60500 	
03700-U			Retinoscope 18100 





















Replacement Lamps

Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
03800-U		<p>PanOptic Ophthalmoscope</p>  <p>11810 11820</p>	<p>PanOptic Ophthalmoscope</p>  <p>11800</p>
03900-U		<p>Ophthalmoscopes</p>  <p>13000</p>	<p>Ophthalmoscopes</p>  <p>12810 12810</p>
04100-U			<p>Examination Lights</p>  <p>48400 48410</p>
04200-U	 <p>(shown at 25% of actual size)</p>		<p>Examination Lights</p>  <p>48600 48610 48625 48635 48700</p>
04400-U			<p>Ophthalmoscopes</p>  <p>11400 11411 11470/11475</p>
04430-U	<p>LS135</p>  <p>44300 44310</p>		
04450-U	<p>LS150</p>  <p>44500 44510</p>		
04500-U			<p>Retinoscope</p>  <p>18010</p>
04600-U			<p>Retinoscope</p>  <p>17710</p>







Replacement Lamps

Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
04700-U		<p>Standard Laryngoscopes</p> <p>680 Series Sizes 0-1 </p> <p>690 Series Sizes 1-2 </p> <p>692 Series Sizes 1-2 </p>	<p>Standard Laryngoscopes</p> <p>634 Series Sizes 0-1 </p> <p>684 Series Sizes 0-1 </p> <p>664 Series Size 1 </p> <p>674 Series Size 1 </p> <p>694 Series Sizes 1-2 </p>
04800-U		<p>Standard Laryngoscopes</p> <p>680 Series Sizes 2-4 </p> <p>690 Series Sizes 3-4 </p> <p>692 Series Sizes 3-4 </p>	<p>Standard Laryngoscopes</p> <p>634 Series Sizes 2-4 </p> <p>640 Series Size 1 </p> <p>640 Series Sizes 2-4 </p> <p>664 Series Sizes 2-4 </p> <p>674 Series Sizes 2-4 </p> <p>684 Series Sizes 2-4 </p> <p>694 Series Sizes 3-4 </p> <p>Transilluminator 40415 </p>
04900-U		<p>Ophthalmoscopes</p> <p> 11720  11730/11735</p>	<p>Ophthalmoscopes</p> <p> 11620  11630</p>
06000-U		<p>Laryngoscope Handles</p> <p>   </p> <p>60813 60814 60713 60815</p>	<p>Laryngoscope Handles</p> <p> </p> <p>60803 60804</p>
06100-U			<p>Examination Lights</p> <p>48400/48435 </p>
06200-U		<p>AudioScope</p> <p>23300 </p>	<p>AudioScope</p> <p>23000, 23020, 23040 </p>
06300-U			<p>Examination Lights</p> <p> 44100  44110</p>
06400-U			<p>Examination Light</p> <p>44200 </p>

Replacement Lamps

Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
06500-U		MacroView Otoscope 23810, 23820 	
07200-U			Retinoscope 17200 
07600-U		PenLite 76600 	
07800-U		Sigmoidoscope, Anoscope Light Handle 73211 	Sigmoidoscope, Anoscope Light Handle 73210 Vaginal Specula Illuminator 78000 Sigmoidoscope, Anoscope Light Handle 73200   
08200-U		Streak Retinoscope 18245 	Streak Retinoscopes 18200 18235  
08300-U		Spot Retinoscopes 18300 	
08400-U		Retinoscope 93300	
08500-U		LumiView 20520 	




















Replacement Lamps

Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
08800-U		 78800	
09500-U	 (shown at 25% of actual size)	 49500	
09800-U		 88000A 89001A 89000A 88007 88001A	

BATTERIES

REPLACEMENT PARTS

Replacement Batteries

Battery Number	Replacement Batteries	Current Instruments	Absolute Instruments
71960		Lithium Ion Handle 71900 	
72000		70700, 71500 	60700 
72200		71670 	617 Handles 717 Handles 
72240		Binocular Indirect Ophthalmoscope 12000 	71000, 71020 
72300		71000-A 71000-C 	
72420		SureSight 14010 	
72600		PocketScope 72800 	
72610		CompacSet 93400 	

BLOOD PRESSURE CUFF REFERENCE GUIDE

REPLACEMENT PARTS

Manual Device Cuff Part Numbers



	Wall, Desk, Mobile & Pocket Aneroids					Handheld Aneroids	
Connector	1 Tri-Purpose & One-Tube for Bulb					1 Tri-Purpose	
							
Type	Durable		Disposable			Durable	Disposable
Material	Nylon		Soft		Vinyl	Nylon	Soft
Bulb	Without Bulb	With Bulb	Without Bulb	With Bulb	With Bulb	Without Bulb	Without Bulb
Unit of Measure	Each	Each	20 Pack	20 Pack	20 Pack	Each	20 Pack
Sm. Infant	N/A	N/A	N/A	N/A	N/A	Reuse-06-1TP	Soft-06-1TP
Infant	Reuse-07-2TP	Reuse-07-2BV	Soft-07-2TP	Soft-07-2BV	Vinyl-07-2BV	Reuse-07-1TP	Soft-07-1TP
Sm. Child	Reuse-08-2TP	Reuse-08-2BV	Soft-08-2TP	Soft-08-2BV	Vinyl-08-2BV	Reuse-08-1TP	Soft-08-1TP
Child	Reuse-09-2TP	Reuse-09-2BV	Soft-09-2TP	Soft-09-2BV	Vinyl-09-2BV	Reuse-09-1TP	Soft-09-1TP
Sm. Adult	Reuse-10-2TP	Reuse-10-2BV	Soft-10-2TP	Soft-10-2BV	Vinyl-10-2BV	Reuse-10-1TP	Soft-10-1TP
Adult	Reuse-11-2TP	Reuse-11-2BV	Soft-11-2TP	Soft-11-2BV	Vinyl-11-2BV	Reuse-11-1TP	Soft-11-1TP
Adult Long	Reuse-11L-2TP	Reuse-11L-2BV	Soft-11L-2TP	Soft-11L-2BV	Vinyl-11L-2BV	Reuse-11L-1TP	Soft-11L-1TP
Lg. Adult	Reuse-12-2TP	Reuse-12-2BV	Soft-12-2TP	Soft-12-2BV	Vinyl-12-2BV	Reuse-12-1TP	Soft-12-1TP
Lg. Adult Long	Reuse-12L-2TP	Reuse-12L-2BV	Soft-12L-2TP	Soft-12L-2BV	Vinyl-12L-2BV	Reuse-12L-1TP	Soft-12L-1TP
Thigh	Reuse-13-2TP	Reuse-13-2BV	Soft-13-2TP	Soft-13-2BV	Vinyl-13-2BV	Reuse-13-1TP	Soft-13-1TP

	DS44		DS45	
Connector	Connector Integrated Into Cuff			
				
Type	Durable		Reusable	
Material	Nylon		Nylon	
Bulb	Without Bulb	With Bulb	Without Bulb	With Bulb
Unit of Measure	Each	Each	Each	Each
Child	DS44-DUR09	DS44-DUR09F	DS45-DUR09	DS45-DUR09F
Child Print	DS44-DUR09P	DS44-DUR09PF	DS45-DUR09P	DS45-DUR09PF
Sm. Adult	DS44-DUR10	DS44-DUR10F	DS45-DUR10	DS45-DUR10F
Adult	DS44-DUR11	DS44-DUR11F	DS45-DUR11	DS45-DUR11F
Lg. Adult	DS44-DUR12	DS44-DUR12F	DS45-DUR12	DS45-DUR12F

BLOOD PRESSURE CUFF REFERENCE GUIDE


REPLACEMENT PARTS

Electronic Device Cuff Part Numbers

	Atlas, Vital Signs and Spot Vital Signs Devices (except LXi)			Propaq and PIC Monitors		
Connector	1 Screw-Type			1 Locking Type		
						
Type	Reusable	Disposable		Reusable	Disposable	
Material	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl
Unit of Measure	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack
Sm. Infant	Reuse-06-1SC	Soft-06-1SC	Vinyl-06-1SC	Reuse-06-1MQ	Soft-06-1MQ	Vinyl-06-1MQ
Infant	Reuse-07-1SC	Soft-07-1SC	Vinyl-07-1SC	Reuse-07-1MQ	Soft-07-1MQ	Vinyl-07-1MQ
Sm. Child	Reuse-08-1SC	Soft-08-1SC	Vinyl-08-1SC	Reuse-08-1MQ	Soft-08-1MQ	Vinyl-08-1MQ
Child	Reuse-09-1SC	Soft-09-1SC	Vinyl-09-1SC	Reuse-09-1MQ	Soft-09-1MQ	Vinyl-09-1MQ
Sm. Adult	Reuse-10-1SC	Soft-10-1SC	Vinyl-10-1SC	Reuse-10-1MQ	Soft-10-1MQ	Vinyl-10-1MQ
Adult	Reuse-11-1SC	Soft-11-1SC	Vinyl-11-1SC	Reuse-11-1MQ	Soft-11-1MQ	Vinyl-11-1MQ
Adult Long	Reuse-11L-1SC	5082-11L-1SC	Vinyl-11L-1SC	Reuse-11L-1MQ	Soft-11L-1MQ	Vinyl-11L-1MQ
Lg. Adult	Reuse-12-1SC	5082-12-1SC	Vinyl-12-1SC	Reuse-12-1MQ	Soft-12-1MQ	Vinyl-12-1MQ
Lg. Adult Long	Reuse-12L-1SC	5082-12L-1SC	Vinyl-12L-1SC	Reuse-12L-1MQ	Soft-12L-1MQ	Vinyl-12L-1MQ
Thigh	Reuse-13-1SC	5082-13-1SC	Vinyl-13-1SC	Reuse-13-1MQ	Soft-13-1MQ	Vinyl-13-1MQ

	Hewlett Packard, Spacelabs and Datascope Monitors			Critikon DiN/Amap Monitors			Spot LXi and GE Marquette Monitors		
Connector	1 Bayonet			2 Screw-Type			2 Locking-Type		
Type	Reusable	Disposable		Reusable	Disposable		Reusable	Disposable	
Material	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl
Unit of Measure	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack
Sm. Infant	Reuse-06-1HP	Soft-06-1HP	Vinyl-06-1HP	N/A	N/A	N/A	N/A	N/A	N/A
Infant	Reuse-07-1HP	Soft-07-1HP	Vinyl-07-1HP	Reuse-07-2SC	Soft-07-2SC	Vinyl-07-2SC	Reuse-07-2MQ	Soft-07-2MQ	Vinyl-07-2MQ
Sm. Child	Reuse-08-1HP	Soft-08-1HP	Vinyl-08-1HP	Reuse-08-2SC	Soft-08-2SC	Vinyl-08-2SC	Reuse-08-2MQ	Soft-08-2MQ	Vinyl-08-2MQ
Child	Reuse-09-1HP	Soft-09-1HP	Vinyl-09-1HP	Reuse-09-2SC	Soft-09-2SC	Vinyl-09-2SC	Reuse-09-2MQ	Soft-09-2MQ	Vinyl-09-2MQ
Sm. Adult	Reuse-10-1HP	Soft-10-1HP	Vinyl-10-1HP	Reuse-10-2SC	Soft-10-2SC	Vinyl-10-2SC	Reuse-10-2MQ	Soft-10-2MQ	Vinyl-10-2MQ
Adult	Reuse-11-1HP	Soft-11-1HP	Vinyl-11-1HP	Reuse-11-2SC	Soft-11-2SC	Vinyl-11-2SC	Reuse-11-2MQ	Soft-11-2MQ	Vinyl-11-2MQ
Adult Long	Reuse-11L-1HP	Soft-11L-1HP	Vinyl-11L-1HP	Reuse-11L-2SC	Soft-11L-2SC	Vinyl-11L-2SC	Reuse-11L-2MQ	Soft-11L-2MQ	Vinyl-11L-2MQ
Lg. Adult	Reuse-12-1HP	Soft-12-1HP	Vinyl-12-1HP	Reuse-12-2SC	Soft-12-2SC	Vinyl-12-2SC	Reuse-12-2MQ	Soft-12-2MQ	Vinyl-12-2MQ
Ld. Adult Long	Reuse-12L-1HP	Soft-12L-1HP	Vinyl-12L-1HP	Reuse-12L-2SC	Soft-12L-2SC	Vinyl-12L-2SC	Reuse-12L-2MQ	Soft-12L-2MQ	Vinyl-12L-2MQ
Thigh	Reuse-13-1HP	Soft-13-1HP	Vinyl-13-1HP	Reuse-13-2SC	Soft-13-2SC	Vinyl-13-2SC	Reuse-13-2MQ	Soft-13-2MQ	Vinyl-13-2MQ

Electronic Device Cuff Part Numbers















Spot Vital Signs LXi and GE Marquette Monitors			
Connector	2 Locking-Type		
			
Type	Reusable	Disposable	
Material	Nylon	Soft	Vinyl
Unit of Measure	Each	20 Pack	20 Pack
Sm. Infant	N/A	N/A	N/A
Infant	Reuse-07-2MQ	Soft-07-2MQ	Vinyl-07-2MQ
Sm. Child	Reuse-08-2MQ	Soft-08-2MQ	Vinyl-08-2MQ
Child	Reuse-09-2MQ	Soft-09-2MQ	Vinyl-09-2MQ
Sm. Adult	Reuse-10-2MQ	Soft-10-2MQ	Vinyl-10-2MQ
Adult	Reuse-11-2MQ	Soft-11-2MQ	Vinyl-11-2MQ
Adult Long	5082-11L-2MQ	Soft-11L-2MQ	Vinyl-11L-2MQ
Lg. Adult	5082-12-2MQ	Soft-12-2MQ	Vinyl-12-2MQ
Ld. Adult Long	5082-12L-2MQ	Soft-12L-2MQ	Vinyl-12L-2MQ
Thigh	5082-13-2MQ	Soft-13-2MQ	Vinyl-13-2MQ

Neonatal Monitors		
Connector	Luer Slip	
Tubes	One-Tube	Two-Tube
Type	Disposable	
Material	Soft	
Unit of Measure	10 Pack	10 Pack
Neo 1	5082-101-1	5082-101-2
Neo 2	5082-102-1	5082-102-2
Neo 3	5082-103-1	5082-103-2
Neo 4	5082-104-1	5082-104-2
Neo 5	5082-105-1	5082-105-2






BLOOD PRESSURE CONNECTORS

REPLACEMENT PARTS

Replacement Connectors

MALE		FEMALE	
Screw-Type			
5082-164 	Plastic Male Screw-Type with Female Luer Slip Inside, 5/32" (4 mm) Hose Barb. For Welch Allyn Atlas and Vital Signs Monitors, Welch Allyn Spot, and Critikon cuffs, 10/pack	5082-161 	Metal Female Screw-Type, 1/8" (3.2 mm) Hose Barb. For Welch Allyn Atlas and Vital Signs Monitors, Welch Allyn Spot, and Critikon monitors, 10/pack
		5082-174 	Plastic Female Screw-Type, 5/32" (4 mm) Hose Barb. For Welch Allyn Atlas and Vital Signs Monitors, Welch Allyn Spot, and Critikon monitors, 10/pack
		5082-175 	Plastic Female Screw-Type, 1/8" (3.2 mm) Hose Barb. For Welch Allyn Atlas and Vital Signs Monitors, Welch Allyn Spot, and Critikon monitors, 10/pack
Luer Locks			
5082-165 	Plastic Male Luer Lock, 5/32" (4 mm) Hose Barb, 10/pack	5082-168 	Plastic Female Luer Lock, 5/32" (4 mm) Hose Barb. For Sphygmomanometer cuffs; BCI, Colin, and Datascope Passport XG and Expert monitor cuffs, 10/pack
Luer Slips			
5082-166 	Plastic Male Luer Slip, 5/32" (4 mm) Hose Barb, 10/pack	5082-164 	Plastic Male Screw-Type with Female Luer Slip Inside, 5/32" (4 mm) Hose Barb, 10/pack
5082-167 	Metal Male Luer Slip, 5/32" (4 mm) Hose Barb, 10/pack	5082-169 	Metal Female Luer Slip, 5/32" (4 mm) Hose Barb, 10/pack
5082-176 	Plastic Male Luer Slip, 3/32" (2.4 mm) Hose Barb. For neonatal cuffs, 10/pack		
Locking-Type			
5082-172 	Plastic Male Locking Adapter, 1/8" (3.2 mm) Hose Barb. For Welch Allyn Propaq, GE Marquette two-tube, and Datex-Ohmeda* monitors, 10/pack	5082-182 	Plastic Female Locking Connector, 5/32" (4 mm) Hose Barb. For Welch Allyn Propaq, GE Marquette two-tube, and Datex-Ohmeda* monitors, 10/pack
5082-178 	Plastic Male Locking Adapter, 5/32" (4 mm) Hose Barb. For Welch Allyn Propaq, GE Marquette two-tube, and Datex-Ohmeda* monitors, 10/pack		
*Datex-Ohmeda uses one male and one female connector on its two-tube monitors and cuffs.			






Replacement Connectors

MALE		FEMALE	
Bayonets			
5082-179 	Metal Male Bayonet, 5/32" (4 mm) Hose Barb. For HP/ Philips, Siemens, Spacelabs, Criticare, IVAC/Alaris 4415, and Datascope Accutorr Plus and Passport 2 monitors, 10/pack	5082-177 	Metal Female Bayonet, 1/8" (3.2 mm) Hose Barb. For HP/ Philips, Siemens, Spacelabs, Criticare, IVAC/Alaris 4415, and Datascope Accutorr Plus and Passport 2 monitors, 2/pack
5082-184 	Plastic Male Bayonet, 5/32" (4 mm) Hose Barb. For HP/ Philips, Siemens, Spacelabs, Criticare, IVAC/Alaris 4415, and Datascope Accutorr Plus and Passport 2 monitors, 10/pack	5082-185 	Metal Female Bayonet, 5/32" (4 mm) Hose Barb. For HP/ Philips, Siemens, Spacelabs, Criticare, IVAC/Alaris 4415, and Datascope Accutorr Plus and Passport 2 monitors, 2/pack
Tri-Purpose Sphygmomanometer Connector			
5082-200 	Plastic Male M9 thread with Female luer lock and Female luer slip inside, 5/32" (4 mm) Hose Barb, 1 each		








TUBING AND ACCESSORIES

REPLACEMENT PARTS

Replacement Tubes and Accessories

PART NUMBER	DESCRIPTION
Straight Tubing with Connectors and Adapters	
<p>5082-153</p> 	<p>Sphyg tubing with 2 plastic male luer slip cuff connectors on one end and plastic hose barb for sphyg tubing and 1/3" (33 cm) tubing for inflation bulb on the other, 10/pack</p>
<p>5082-244</p> 	<p>Monitor tubing, 10 ft (3 m) with metal female bayonet adapter connector, for Welch Allyn Propaq and HP/Philips neonatal monitors, 1/pack</p>
<p>5082-245</p> 	<p>Monitor tubing, 10 ft (3 m) with metal female luer slip adapter. For Welch Allyn Propaq and HP/Philips monitors, 1/pack</p>
<p>5082-246</p> 	<p>Monitor tubing, 10 ft (3 m) with metal female screw-type adapter, for Welch Allyn Propaq and HP/Philips monitors, 1/pack</p>
<p>5082-247</p> 	<p>Monitor tubing, 9 ft 2" (2.8 m) with plastic male screw-type connector (female luer slip inside), for Spacelabs neonatal monitors, 1/pack</p>

Replacement Tubes and Accessories

PART NUMBER		DESCRIPTION
Y-Tubing with and without Connectors and Adapters		
5082-148		Y-tubing with two plastic female screw-type adapters on one end and one plastic bayonet connector on the other, 10/pack
5082-149		Y-tubing with two plastic female screw-type adapters on one end and an 1/8" (3.2 mm) I.D. open end on the other, 10/pack
5082-154		Y-tubing with two plastic female screw-type connectors on one end and one metal female bayonet adapter on the other, 2/pack
5082-183		Y-tubing without connectors. 1/8" (3.2 mm) I.D., 10/pack
Tube Connectors		
5082-170		Plastic connector for 5/32" (4 mm) tubing, 10/pack
5082-171		Plastic connector for 1/8" (3.2 mm) to 5/32" (4 mm) I.D. tubing, 10/pack
Other Accessories		
5082-163		Plastic washing plug (closed end) for 5/32" (4 mm) I.D. tubing, 10/pack