

Welch Allyn 1500 Patient Monitor Competency Checklist

Name (print): _____ Unit: _____ Date: _____

- _____ Turn the Welch Allyn 1500 Monitor on; describe the steps to turn the device off.
- _____ Obtain a manual NIBP. Change the NIBP to automatic mode; then change the time intervals.
- _____ Demonstrate how to silence a patient alarm. What are the alarms silence time options?
- _____ Describe how to change the arrhythmia settings at the device.
- _____ Describe the steps to obtain a 12-Lead Resting ECG.
- _____ Explain the process of zeroing the IBP. Which pressure connection determines the HR/PR?
- _____ Display patient vital sign trends; how many hours of trends can the monitor store?
- _____ Describe the Patient Standby Mode. What happens when you resume patient monitoring?
- _____ (True or False) Depending on the temperature sensor type, the sensor can be applied to the ear, the skin, or to the rectum.
- _____ (True or False) The monitor can be cleaned using the following: 70% solution isopropyl alcohol, neutral mild detergent solution, all products designed for cleaning plastic.