

Braun ThermoScan® PRO 4000 Competency Checklist

Name (print): _____ Unit: _____ Date: _____

- ___ Demonstrate how to place a probe cover on the thermometer.
- ___ True or False. Always use a new, clean, Welch Allyn probe cover because used ones will result in inaccurate readings.
- ___ Review the steps to clean a dirty thermometer lens.
- ___ Demonstrate proper placement of the probe into the ear canal and take an ear temperature.
- ___ True or False. A long beep and a solid green light signal the end of a measurement with stable probe placement in the ear.
- ___ Press the 'I/O' button to recall the last temperature taken. What was the reading?
- ___ List four factors that MAY interfere with accurate tympanic temperatures.
- ___ Demonstrate how to replace the 2 AA batteries in the thermometer.