



Vital Signs Monitor 300 Series | Automated Vital Signs Monitor

Welch Allyn Vital Signs Monitor 300 Series

Advanced Monitoring Made Affordable



WelchAllyn[®]

Advancing Frontline Care[™]

Welch Allyn Vital Signs Monitor 300 Series.

Simple. Flexible. Affordable.

Your spot-check and continuous vital signs monitor, all in one.

As a clinician you are under pressure to improve outcomes, boost operating efficiencies and lower costs every day. It's a daunting challenge, but one that can be made a little easier thanks to the Welch Allyn Vital Signs Monitor 300 Series – an affordable vital signs monitor that's both flexible and easy to use.

Like the original Welch Allyn Vital Signs Monitor, the 300 Series helps enhance your productivity by providing automated vital signs before, during and after procedures. The simple, push-button menu operation and bright display allow you to quickly learn and easily obtain the important information you need, so you spend more time with your patients and less time with your equipment.

In addition to the great features you've come to expect from Welch Allyn, the Vital Signs Monitor 300 Series provides these exciting features:

- Separate patient modes with individual settings for adult, paediatric and neonatal support
- Fast, accurate temperature readings with industry-leading Welch Allyn SureTemp® Plus thermometry
- Simple, straightforward operation and intuitive menu navigation with our innovative monitor design

And because the 300 Series monitor is ideal for both routine assessment and continuous monitoring of vital signs, it's like getting two devices at one low price.



Procedures

The Vital Signs Monitor 300 Series is the ideal solution for taking routine patient vital signs or monitoring before, during and after procedures, minor surgery and treatments in a variety of settings, including:

- Medical/Surgical Areas
- Paediatric Wards
- Emergency Departments
- Hospital Special Procedures Departments
- Labour and Delivery Rooms
- Outpatient Surgery Departments
- Outpatient Clinics
- Specialty Physician Clinics
- Urgent Care Units
- Dialysis Units
- Dental Surgeries

Comprehensive Features

Non-invasive Blood Pressure

The Welch Allyn Vital Signs Monitor 300 Series begins as a non-invasive blood pressure (NIBP) monitor. Even in its most basic configuration, it delivers the features you need.

- Multiple Patient Modes – Adult, Paediatric and Neonatal operating modes
- Automatic Blood Pressure Mode – choose from 13 measurement intervals
- Programmable Alarms – for high/low systolic, diastolic and pulse rate
- LCD Display – clock window and menu navigation
- LED Display – bright, clear viewing of vital signs readings
- Built-In Memory – save up to 99 sets of measurement data

SureTemp® Plus Thermometry

Add SureTemp Plus Electronic Thermometry and take advantage of the revolutionary measurement capability developed by Welch Allyn. Obtain oral, rectal and axillary temperatures with the speed you want and the accuracy you demand.

- Oral – approx. 4-6 seconds
- Axillary – Adults, approx. 15 seconds; paediatrics, approx. 10 seconds
- Rectal – approx. 10 seconds



Pulse Oximetry

The Vital Signs Monitor 300 Series can be configured with motion-tolerant pulse oximetry (SpO₂) technology from Masimo® or *OxiMax*® technology from Nellcor.


- Programmable alarms – high/low saturation levels
- Wide variety of patient sensors
- Optional pulse tone feature – audible heart rate and SpO₂ saturation indicator

Optional Printer

Add an integrated thermal printer for trend analysis or record keeping and get what you need at the touch of a button.

Simple and Straightforward

The Welch Allyn Vital Signs Monitor 300 Series has a very short learning curve. Simple push-button programming and a large, digital display make learning and using the monitor easy and hassle-free.



Take It Wherever You Need It

Carry it, set it on a countertop, wheel it on its own dedicated mobile stand, or mount it on a wall. Do it all – switching from one flexible mounting position to another is simple. The Welch Allyn Vital Signs Monitor 300 Series goes wherever your patient goes – effortlessly.

Ordering Information

53000	NIBP only (includes Pulse Rate and MAP)
5300P	NIBP/Printer
530T0	NIBP/Temperature
530TP	NIBP/Temperature/Printer
53S00	NIBP/Masimo Pulse Oximetry
53S0P	NIBP/Masimo Pulse Oximetry/Printer
53ST0	NIBP/Masimo Pulse Oximetry/Temperature
53STP	NIBP/Masimo Pulse Oximetry/Temperature/Printer
53N00	NIBP/Nellcor Pulse Oximetry
53N0P	NIBP/Nellcor Pulse Oximetry/Printer
53NT0	NIBP/Nellcor Pulse Oximetry/Temperature
53NTP	NIBP/Nellcor Pulse Oximetry/Temperature/Printer
008-0835-01	Mobile Stand with Accessory Basket
008-0834-01	Wall Mount with Accessory Basket
008-0890-00	IV Pole Mount

Language version suffix:

E=English, F=French, G=German, I=Italian, S=Spanish,
C=Chinese, D=Dutch, P=Polish, U=Portuguese, V=Swedish

Transformer version suffix:

1=North American, 2=European, 4=UK, 6=Australian

UK Customer Services

Welch Allyn UK Ltd, Cublington Road, Aston Abbots, Buckinghamshire HP22 4ND, England

Tel: +44 (0) 207 365 6780 Fax: +44 (0) 207 365 9694

South Africa Customer Services

Welch Allyn South Africa (PTY) Ltd, 394 Surrey Avenue, Ferndale, Randburg 2194, South Africa

P O Box 809, Fontainebleau 2032, South Africa

Tel: +27 (0) 11 777-7555 Fax: +27 (0) 11 777-7556

Welch Allyn European Customer Services,

Navan Business Park, Dublin Road, Navan, Co. Meath, Republic of Ireland

Tel: +353 46 90 67790 Fax: +353 46 90 67754

www.welchallyn.com

©2005 Welch Allyn. All rights reserved. SureTemp is a registered trademark of Welch Allyn, Inc.
Nellcor and OxiMax are registered trademarks of Nellcor Puritan Bennett Inc. Masimo is a trademark
of Masimo Corporation.

SM2615EU RevD



WelchAllyn®

Advancing Frontline Care™